






*Your medicine dry, clean and close  
at hand*




*The tablets or capsules are stored in  
a plastic inner tube*



*Ridges on the lid make the medi-  
cation accessible with a single  
slide of the thumb*



*Holds up to 10 smaller tablets or two  
average-sized capsules*



*The Ingrid products are made in  
Sweden*

## INGRID Med Tablet Tube

Ingrid Med Tablet Tube safely stores medications, making them readily accessible – wherever you go.

Simply attach the Tablet Tube to your key ring or slip it in your pocket or purse

### Plastic inner tube and aluminum outer shell

The tablets or capsules are stored in a plastic inner tube, which is encased, in a sturdy, scratch-resistant and waterproof aluminum outer shell.

The plastic tube has two components – a canister which houses the tablets or capsules and a lid. Ridges on the lid make the medication accessible with a single slide of the thumb.

Never store medications directly in the aluminum shell without using the plastic inner tube.

### 10 smaller tablets or two average-sized capsules

Ingrid Med Tablet Tube holds up to 10 smaller tablets (up to 8mm Ø) or two average-sized capsules (up to 18mm long / 8mm Ø).

### Handling

To properly preserve your medications, and keep moisture out, always fully tighten the lid. Do not expose the Tablet Tube to high temperatures (such as a locked car). So as not to risk expiration, regularly change the stored medications. Please consult your doctor or pharmacist if you are uncertain whether to store certain medications in the Ingrid Med Tablet Tube.

Ingrid Med Tablet Tube is not intended to store medications used to treat acute, life-threatening conditions.